

CRISPR knock-in protocol for cell lines (Ayaz Najafov)

(based on the 1% DMSO knock-in efficiency enhancement discovery by Stratigopoulos et al., PLoS One. 2018 Jun 4;13(6):e0198637.)

- 1) Seed HT-29 cells into 6-well plates in 3ml McCoy's medium (10% FBS + Antibiotic/Antimycotic).
- 2) 24 hrs later, at 80% confluence, treat with 1% DMSO (in the old medium).
- 3) 24 hrs after DMSO addition transfect the cells (in the old medium).
 - Per well:**
 - **Tube A:** 125 µl OptiMEM + 2 µg pX459-puro-sgRNA plasmid + 3 µg Mutant PCR donor + 5 µl P3000
 - **Tube B:** 125 µl OptiMEM + 7.5 µl Lipofectamine 3000
 - Mix the tubes A and B. Pipette to mix. Incubate at 25°C for 15min. Deliver 250 µl per well.
- 4) 24 hrs later (1 day post-transfection), expand each well into a 15cm dish.
- 5) 24 hrs later (2 days post-transfection), add fresh medium with puromycin (2.5 µg/ml).
- 6) 48 hrs later (4 days post-transfection), add fresh medium with puromycin (2.5 µg/ml).
- 7) Single-cell sort by flow cytometry, do limited dilution or keep in 15cm plates for colony picking.
- 8) At 7 days post-transfection, switch to puromycin-free medium.
- 9) Grow for 2 weeks.
- 10) Make two replicates of the clones, in 12-well plates (one for propagation and one for genotyping). Spin the plates to enhance cell attachment: 1,000xg, 3min, 25°C.
- 11) 48 hrs later:
 - a. Feed with the replicates for propagation with 1ml new medium
 - b. Trypsinize the replicates for genotyping and transfer the cell pellets into 8-strip PCR tubes.
- 12) Generate cell lysates using Tail Lysis protocol (Truett et al., BioTechniques 29:52-54 (July 2000))
 - a. Add 75 µl of tail lysis buffer (25 mM NaOH, 1 mM EDTA).
 - b. Heat at 95°C for 1hr.
 - c. Put tubes on ice for 5min.
 - d. Add 75 µl of neutralization buffer (40 mM Tris-HCl, pH5, dissolve Tris-HCl, no pH adjustment).
- 13) Setup PCR:
 - a. 5 µl 2X EmeraldAmp® GT PCR Master Mix
 - b. 5 µl Nuclease-free water
 - c. 0.5 µl 10 µM F+R primer stock
 - d. 0.5 µl of the cell lysate.
 - 98°C 2min (hot start)
 - 40 cycles:
 - 98°C 15sec
 - 60°C 30sec
 - 72°C 1min
 - 72°C 5min
- 14) Prepare the PCR products for sequencing:
 - a. 11 µl PCR product + 1 µl Exol + 2 µl rSAP
 - b. Incubate at 37°C for 30min.
 - c. Incubate at 80°C for 20min.
 - d. Cool to 4°C.
- 15) Send the whole reaction for sequencing with the forward primer used to generate the PCR donor.

Materials:

EmeraldAmp® GT PCR Master Mix	Clontech, RR310B
Exonuclease I (Exol)	NEB, M0293S
Shrimp Alkaline Phosphatase (rSAP)	NEB, M0371S

Mutant PCR donor = ~150bp PCR amplicon (with the desired point mutation in the middle) purified using Qiagen's PCR purification kit.