

# Purification of GST-Fusion Protein

## Materials:

<b>PBS:</b>	140 mM NaCl, 2.7 mM KCl, 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 1.8 mM KH <sub>2</sub> PO <sub>4</sub> (pH 7.3)
<b>PBST:</b>	PBS + 0.5% Triton X-100
<b>Glutathione Resin:</b>	Glutathione Sepharose 4B (Pharmacia 17-07556-01)
<b>Elution Buffer:</b>	100 mM Tris (pH 8.5-9.0), 20 mM reduced glutathione (aliquot and store @ -20 deg C)
<b>0.5 M IPTG</b>	prepare in H <sub>2</sub> O and freeze in aliquots
<b>Protease Inhibitors</b>	add to solutions immediately before use

## Procedure:

1. Induce expression of GST-fusion protein
  - ✦ grow 5 ml O/N culture of E. coli strain harboring pGEX plasmid in LB + Carb (50 ug/ml)
  - ✦ inoculate 500 mls LB+Carb media with 5 mls saturated O/N culture
  - ✦ grow culture at 30 degrees until mid-log phase (OD<sub>600</sub> ~ 0.6 - 0.7)
    - ✦ typically 2.5-3.5 hours for TG1 host strains
    - ✦ collect 100 ul cells for analysis and resuspend in 10 ul PBS
  - ✦ induce fusion protein expression by adding IPTG to 0.1 mM
    - ✦ grow cells for 3 additional hours at 30 degrees (final OD<sub>600</sub> ~ 1.2 - 1.3)
    - ✦ collect 100 ul cells for analysis and resuspend in 20 ul PBS buffer
  - ✦ harvest cells by centrifugation 3000 x g<sub>av</sub> for 8 min (5K rpm in SLA-1500 rotor)
  - ✦ resuspend cells in 20 mls PBS
  - ✦ concentrate cells by centrifugation 3000 x g for 8 min and resuspend in final volume of PBS no greater than 5 mls
  - ✦ transfer to 50 ml conical tube
2. Prepare extract (all work done @ 4 degrees C)
  - ✦ freeze cells at -80 degrees 1 hr - O/N
  - ✦ thaw cells and add protease inhibitor cocktail
  - ✦ lyse cells by sonication: 4 - 6 x 10 seconds with 45 second rests on ice between bursts
    - ✦ (Fisher 550 Sonic Dismembrator @ setting 4.5 with microtip)
    - ✦ monitor cell lysis by examination with compound microscope
  - ✦ add Triton X-100 to final concentration of 1% and rock gently for 20 minutes
  - ✦ clarify extract by centrifugation
    - ✦ transfer lysate to centrifuge tube
    - ✦ wash and collect residual with 5 mls GST buffer and combine with lysate
    - ✦ spin down insoluble material 10 min @ 8000 x g<sub>av</sub> (10K rpm in SS-34 rotor)
  - ✦ transfer supernatant to clean tube
    - ✦ save 10 ul for analysis
  - ✦ resuspend pellet in 10 mls PBST and save 10 ul for analysis
3. Prepare glutathione resin
  - ✦ add 1.33 mls 75% glutathione sepharose slurry to 15 ml conical tube (= 1 ml bed volume)
  - ✦ spin down resin 2 min @ 500 x g (2K rpm in IEC tabletop centrifuge) and remove supernatant
  - ✦ wash 2x with cold PBS

- ⌘ wash 1x with cold PBST
- 4. Bind GST-fusion protein to GSH-resin
  - ⌘ transfer clarified extract to prepared glutathione resin
  - ⌘ rock gently @ 4 degrees for 2 hours to allow binding
  - ⌘ remove unbound material by centrifugation 2 min @ 500 x g and removal of supernatant
  - ⌘ wash resin 2x with cold PBST
  - ⌘ wash resin 1x with cold PBS
- 5. Elute bound protein from resin
  - ⌘ transfer resin to empty column and allow resin to settle
  - ⌘ elute protein from resin in 10 x 250 ul fractions
  - ⌘ determine fractions containing protein by Bradford assay
    - ⌘ GST-fusion proteins typically eluted in fractions 3-6
  - ⌘ asses protein quality by SDS-PAGE of protein containing fractions
  - ⌘ pool protein containing fractions and store at -80 degrees

### **Notes & Misc:**

- ⌘ Often a 70 kD protein co-purifies with the GST-fusion protein. This is likely the chaparonin DnaK and can sometimes be removed by treating the clarified lysate with 2 mM ATP, 10 mM MgSO<sub>4</sub> before binding the GST-fusion to the resin.
- ⌘ Use 25-30 degrees instead of 37 degrees for cultures to improve solubility of GST-fusion proteins
- ⌘ 1 ml bed volume of resin is reported to bind approx 5-8 mg protein

### **References:**

- ⌘ GST Gene Fusion System (Pharmacia Biotech)
- ⌘ solubilization from inclusion bodies: Frangioni, J.V. and B.G. Neel (1993). Anal. Biochem. 210:179-187.

*protocol compiled by Chad Rappleye*

*Aroian Lab Protocols*